

7-1-96

MIDI Prep of Plasmid DNA from Transformed DM1

- 1) culture started from colony w/ 100 μ l cells (see previous) plated @ 3 μ ml \rightarrow then 30 ml
- 2) MIDI Prep of plasmid DNA (see Qiagen MIDI Prep protocol)
- 3) DNA resuspended in \sim 30 μ l TE

Notes: * culturing of DM1 was very slow; about 8 hr. to get from clear media to turbid
* centrifugation to pellet DNA during MIDI prep was not favorable; DNA not readily visible so DNA ~~may~~ (or sufficient quantities of DNA) may not have been isolated


7-1-96 → 7-2-96

Digestion of pEGFP-N w/ Age-I & Bsp EI and pFOXSB w/ Bsp EI

total rxn-
volumes100 μ l7-2-96

- 1) pEGFP-N (³⁰~~25~~ μ g from previous experiment) was digested w/ Age I @ room temp. overnight (4 μ l enzyme)
- 2) DNA was reprecipitated w/ cold EtOH, resuspended and cut w/ Bsp EI (4 μ l) for ~4 hr. @ 37°C
- 3) Concurrently, pFOXSB (5 μ g) was cut w/ Bsp EI (4 μ l) and incubated @ 37° for ~4 hr.
- * pEGFP-N (cut w/ Bsp EI) and pFOXSB was digested in 50 μ l total volume rxn.
- 4) Digested pFOXSB was then mixed w/ CIP for ²⁵~~20~~ minutes

- 5) Samples run on 1.5% agarose gel (100 V) w/ λ marker



* used only
1/2 amt. of
EtBr in
agarose

pFOXSB

pEGFP-N

 λ

cut w/:

Bsp EI

Bsp EI

+ Age I

7-2-96

6) pFOXSB (BspEI cut) band and pEGFP-N 750 bp band was cut out from gel and purified using Qiagen Gel Extraction Kit

* pFOXSB band = 215 mg

* pEGFP-N 750 bp band = 206 mg

7) DNA eluted from spin column was resuspended in 50 μ l 1x TE

Ligation of pFOXSB and pEGFP-N 750 bp fragment (EGFP + NLS)

1) The following ligation reactions were setup:

<u>Tube #</u>	<u>contents</u>
1	1 μ l cut pFOXSB + 5 μ l pEGFP-N fragment + 1 μ l ligase + 1 μ l ligase buffer + 2 μ l ddH ₂ O (ligase added last) tot. vol = 10 μ l
2	1 μ l cut pFOXSB + 10 μ l pEGFP-N fragment + 1 μ l ligase + 1.5 μ l buffer + 1.5 μ l ddH ₂ O tot. vol = 15 μ l
(control) 3	1 μ l cut pFOXSB + 9 μ l ddH ₂ O (for background) tot. vol = 10 μ l

2) The above contents were mixed (except for ligase) in Eppi & centrifuged briefly; then ligase was added

Ligation (cont'd)

7-2-96

3) ligation rxn. incubated overnight
@ 14-15°C

7-3-96

Plasmid MIDI Prep

* pellet bacteria from culture by spinning @ 40° for 10' speed 4

1) MIDI Prep performed on DM1 transformed cells from a second culture (prepared by Maibe) following culture procedures described in Qiagen plasmid MIDI protocol. (see 7-1-96)

* DM1 culture yielded more cells this time (~3X)

2) for MIDI prep., P1 buffer and P2 buffers were added before resuspension ~~and mixing~~ of cells.

3) After ~~15' incubation in ice~~ the addition of P3 buffer and subsequent incubation on ice for 15', the mixture was aliquoted into Eppi's w/ 1 ml in each tube for a total of 12 tubes.

4) The tubes were then centrifuged @ hi speed for 30' ~~in microcentrifuge~~ in microcentrifuge (@ RT (room temp.))

5) Qiagen tip-100 was equilibrated w/ 4 ml QBT buffer

6) Supernatant from each tube was applied to column

7) Column washed 2x 10 ml QC buffer

8) DNA eluted from column w/ QF buffer (5 ml) and collected in Eppi's (aliquoted to 500 µl per tube)

9) DNA was ppt. by adding 350 µl isopropanol to each tube & centrifuging @ hi speed (RT) for 30'


Plasmid MIDI Prep (cont'd)

7-3-96

10) Ethanol (cold & 70%) was used to wash DNA (500 μ l / tube) and redissolved in 75 μ l TE 1X
 * 75 μ l TE was added to tube 1 and this same TE was used to resuspend ppt. DNA in the rest of the 9 tubes; ultimately, all DNA (originally separated in 10 tubes) is resuspended in 1 tube labelled "pEGFP-N (MIDI) 7-3-96"

run @
100V

11) an aliquot of this prep was then run on an agarose (1%) gel:



pEGFP-N1
pGK-neo (marker)

* 10 μ l λ marker

* pEGFP-N1 & pGK-neo = 1 μ l plasmid

7-3-96

Transformation of Competent cells w/ pFOXEGFP-N1 ligation

*bacteria
used is
HB101

- 1) cells thawed on ice (100ul cells/tube), then add all of ligation rxns. to tubes of cells
- 2) a fourth transformation rxn. was performed using BS ΔNot to transform cells (1ul BS ΔNot) for Maïke
- 3) incubate on ice for 30' ; warm 0.9 ml LB @ 42°
- 4) heat shock @ 42° for 45 sec and ice for 2' ; then add ligation rxns. to ^{LB} cells (all of the rxn. tube)
- 5) incubate @ 37° for 1 hr. ; briefly spin down cells (5 sec.)
- 6) plate on amp plates as follows:
 - a) plate 100ul of 5ul insert, 10ul insert and control ^{remove 900ul & resuspend cells in remaining 100ul media}
 - b) streak plate w/ loop 5ul BS ΔNot
 * 1ul BS ΔNot → 10ul total vol. w/ ddH₂O ; then added to 100ul competent cells
- 7) incubate @ 37° ~~for 1 hr.~~ overnight
- 8) picked 24 colonies from plate w/ 10ul insert and grown in culture tubes w/ 3ml LB overnight in shaker

7-5-96

Miniprep of pFOXEGFP-N1 and Endonuclease Digestion

- 1) put 1 ml of liquid culture (from total of 24 cultures) in Eppi's and spin 50 sec.
- 2) pour off supernatant by ~~flicking~~ inverting tube once; this leaves ~50 μ l in tube; resuspend pellet in remaining media by vortexing
- 3) add 300 μ l TENS (cell lysing) ^{mix} and 150 μ l 3M NaOAc pH 5.2 ^{mix}; centrifuged @ hi speed for 3'
- 4) supernatant transferred to new Eppi's and added 900 μ l 100% cold EtOH; mix by inverting tubes; then, centrifuge @ hi speed for 3'
- 5) aspirate off supernatant and wash w/ cold 70% EtOH; ~~gently invert tubes~~ and centrifuge @ hi speed for 5' (do not mix or invert tubes!)
- a) aspirate off supernatant and resuspended in 30 μ l 1x TE

Digestion

w/ EcoRI

- 1) digestion of pFOXEGFP-N1 w/ EcoRI were setup as follows:


master mix prepared on ice	{	5 μ l DNA (pFOXEGFP-N1; and pFOXSB as control)
		0.1 μ l RNase (10 mg/ml)
		0.5 μ l EcoRI
		1 μ l Buffer I
		3.4 μ l ddH ₂ O
		10 μ l total volume

4

Digestion of pFOXEGFP-N1 (cont'd)

7-5-96

- 2) digestion @ 37° was performed for 2 hrs.
- 3) add 2 μ l loading dye (6X) to each tube and load samples on 1.5% agarose gel and run @ 90 V :



MC12456789101112131415161718192021222324M3


4 → 24

- * 2 and 13 = no DNA
- * something wrong w/ gel?
- * DNA does not appear digested well ; digest > 2hr?
- * not well digested b/c used Buffer 1 V
- * control = pFOXSB (1 μ l plasmid)

Redigestion of pFOXEGFP-N1 miniprep

7-8-96

- 1) digestion rxns. were setup as previously described except used EcoRI buffer (1x in sample). and only ~ 0.5 μ l pFOXSB as control (ran out of pFOXSB)
- 2) digested @ 37° for ~2 hrs. 20 min and loaded onto 1.5% gel @ 100 V :



* 175 band
not seen;
low molar
concentration?!

- * loaded 20 μ l 1x marker
- * control = ~0.5 μ l pFOXSB cut w/ EcoRI (2 μ l)
- * all spl loaded were $\geq 12 \mu$ l
- * lanes between 22 & 23 were empty
- * presence of 850 bp band in 24 but no 175 bp band is observed for all spl; this may have resulted from incomplete digestion and dirty prep
- * will make a new prep of colony 24 and redigest

7-8-96

Reculture of colony 24 for miniprep

- 1) put sterile ~~amp~~ LB (3 ml) into 2 culture tubes each (LB + amp)
- 2) flame an inoculating loop, cool and dip into tube w/ colony 24 and dip into fresh tube w/ LB to reculture; repeat a second time and incubate tubes @ 37° overnight until turbid;
additional tubes (for Maize):

ⓐ BS Δ Not

ⓑ pGK-neo

* for these, inoculated a colony from each plate and cultured overnight in 3 ml LB + amp

Miniprep (Qiagen) of colony 24 & BS Δ Not &

7-9-96

pGK-neo

* bacteria
used is
HB101

- ① pellet bacteria by spinning (speed 4) @ 4° C for 10'; discard supernatant
- ② add 300 μl Buffer P1 and resuspend pellet; transfer to Eppi's
- ③ add 300 μl Buffer P2 gently mix and incubate @ RT for 5'
- ④ add chilled Buffer P3 (300 μl); invert 5-6 x and incubate on ice for 10'
- ⑤ centrifuged @ hi speed for 15'; supernatant moved to 20-tip column (equilibrated)
- ⑥ wash 4 x 1 ml Buffer QC, then DNA eluted w/ 800 μl Buffer QF

Miniprep (cont'd)

7-9-96

⑦ DNA ppt. w/ 560 μ l isopropanol (0.7 vol)
and centrifuged @ hi speed for 30'

⑧ DNA washed w/ 1 ml 70% cold EtOH;
centrifuge 15 sec. and redissolved in 20 μ l 1x TE

Digestion of Colony 24

1) Digestion set up as follows: (2 rxns. made)

1 μ l DNA (pFOXEGFP-N1)

1 μ l EcoRI Buffer

1 μ l EcoRI

7 μ l ddH₂O


10 μ l tot. vol.

Digested for 1 hr. @ 37°C

2) loaded onto 1.5% agarose gel:

* 20 μ l marker

* 12 μ l digestion spl



* put in
EtBr
bath
1 μ g/ml
for 10';
wash for
5'

* where's the 850 bp band?

* funky gel


* ~ 1000 bp band seen


* will redigest

* 10 weight band

observed on

ministry print







Model 373A

Version 1.2.0

Sample 12

Dye Terminator (Any Primer)

Lane 12

Signal: G:409 A:331 T:257 C:164

Points 62

Instrument: 907443

MS 2170

Base 1: 620


Wed, Jul 10, 1996 3:00 PM

X: 0 to 6867 Y: 0 to 1200

Spacing: 10.48

Page 2 of 2

50 JNNNAC GNA 486 T GCA AAC AAN C T NNGGGA AA AA A C T CCA 830 T AACTCCA 840 C AANCC NC



Sac I Digest pFOXEGFP-N1

7-9-96

1) the following digest was setup:

10 μ l DNA (pFOXEGFP-N1)

1 μ l BSA (final conc. 2X)

4 μ l Sac I

5 μ l Buffer I

30 μ l ddH₂O


50 μ l total volume

* a second digestion w/ Eco RI was setup as previously described (1 μ l Eco RI & 1 μ l DNA); both digestions performed in PCR machine @ 37° for 4 hr.

7-10-96

Eco RI digest of pFOXEGFP-N1 was run on 1.5% agarose gel:

* colony
24 #1



* 20 μ l λ marker & 12 μ l rxn. spl

* bad digestion; can see a lot of undigested plasmid; also observe the partial digest band ~ > 1000 bp

→ this may have resulted from incomplete digest to make 850 + 175 \Rightarrow 1025 b band

* 850 bp (GFP & NLS) band observed

CIP of pFOXEGFP-N1

7-10-96

- ① CIP rxn. of Sac I digested vector was setup as follows: (added directly to digestion rxn.)

30 μ l 10X CIP Buffer
219 μ l ddH₂O
1 μ l CIP
+ 50 μ l Sac I digest
300 μ l total vol. ; incubate @ 37° for 30'

- ② add 1 μ l additional CIP, then incubate @ 37° for 30'

- ③ stop rxn. by adding :

(39) 40 μ l ddH₂O
40 μ l 10X STE
20 μ l 10% SDS

400 μ l total volume ; heat @ 68° for 15'

- ④ Extract 2X w/ phenol / chloroform in 1:1 (v/v) ratio

* add phenol / chloroform; vortex; centrifuge @ hi speed for 5'

- ⑤ EtOH precipitate DNA w/ 1 ml cold EtOH ; incubate in dry ice / EtOH bath for 15' ; centrifuge 15' ; wash w/ 1 ml 70% EtOH ; resuspend in 50 μ l 1X TE

15

Ligation of pFOXEGFP-N1 w/ human β -globin intron

7-10-96

1) ligation rxns. were setup as follows:

a) 1 μ l plasmid (pFOXEGFP-N1)
2 μ l insert (β -globin intron)
1 μ l ligase ^{175 bp} (400000 u/ml)
1 μ l ligase buffer (10X)
5 μ l ddH₂O
10 μ l tot. vol.

b) 1 μ l plasmid
6 μ l insert
1 μ l ligase
1.5 μ l lig. buffer
5.5 μ l ddH₂O
15 μ l tot. vol.

c) control: 1 μ l plasmid + 9 μ l ddH₂O = 10 μ l tot. vol.

* 175 bp β -globin human intron was not quantified
2) ligation rxns. were incubated @ 16° overnight

Transformation of HB101 cells w/ pFOXEGFP-N1

7-11-96

Ligation

1) cells thawed on ice and added directly to ligation rxns.; rest of transformation performed as previously described; plate 100 μ l rxn. to each plate (LB + amp)

2) plate on LB+amp (100 μ l) of

a) control b) 2 μ l insert c) 6 μ l insert

* incubate overnight (start 12:00 \Rightarrow 4pm)

7-11-96

7-12-96

Miniprep of pFOXSB

7-11-96

- 1) as described in Sbera's protocol, except put on ice for ~~5'~~ 5' after adding NaOAc (3M pH 5.2)
- 2) after adding 100% EtOH, put in methanol/CO₂ bath for 15'; centrifuge @ hi speed for 3'
- 3) DNA resuspended in 30 μ l 1x TE

Neo I Digestion of primary plasmids used to construct GFP plasmid

* sequence data of pFOXEGFP-N1 showed that the wrong plasmid may have been used to construct GFP plasmid; cutting w/ NcoI should determine whether the right plasmid (pEGFP-C1) and not the wrong old plasmid (pSG5T-C1) was used to construct GFP plasmid. The wrong plasmid pSG5T-C1 will yield an extra low weight fragment (167 bp).

1) Digestion rxn. setup as follows:

Mix made {
5 μ l DNA
0.1 μ l RNase
1 μ l NcoI
1 μ l Buffer 4
2.9 μ l ddH₂O
10 μ l tot. vol.

* for commercial plasmids
pEGFP-C1 and pSG5T-C1
only 1 μ l DNA was used
and added to 4 μ l ddH₂O;
and this was added to 5 μ l
of mix w/ a tot. rxn. vol. of
10 μ l


* also digested minis (pEGFP-N) # 37-48 (positives)
@ 37° for 1 1/2 hrs.; then stored in freezer

* digestion of pEGFP-C1 (prep); pEGFP-N (MIDI prep) and pEGFP-N (unmethylated) was performed

~~Not~~ NcoI Digestion (cont'd)

7-12-96

i) digestion rxns. were run on 2% agarose gel : (120V)



same order:

- | | | |
|--------------------------------------|--------------------------|-------------------|
| 1) Marker | 2) pEGFP-C1 (commercial) | 3) pS65T-C1 (com) |
| 4) pEGFP-C1 (prep) | 5) pEGFP-N (methylated) | |
| 6) pEGFP-N (MIDI prep, unmethylated) | | |

* used 20ul marker

* ~ 12 μ l spl loaded

→ presence of extra 167 bp fragment indicates minis made from old parent vector pSG5T-C1

→ start over w/ pEGFP-C1 prep

7-12-96

Annealing of Oligos

M6528 &
M6529 @
conc. of
1 $\mu\text{g}/\mu\text{l}$

1) added 1 μg (1 μl) of oligo M6528 and M6529 into rxn. w/ 2.5 μl Buffer H (Boehringer-Mannheim) to a total rxn. vol. of 25 μl

2) boiled in water bath for ~1 minute then kept in water bath as it eqs temp. equilibrated to RT

* no kinase

3) put @ 4°C for 2 hr. before using
Digestion of pEGFP-C1 (prep)

1) digestion rxn. setup as follows:

10 μl pEGFP-C1 (prep)

2.5 μl Bsp E1

2.5 μl Bgl II

5 μl Buffer 3

30 μl ddH₂O

50 μl total volume

2) incubate @ 37° for 4 hr.

3) phenol / chloroform and ethanol ppt. the DNA ;
resuspend in 30 μl 1xTE

Ligation of pEGFP-C1 (cut) w/ oligos

7-12-96

MG528 / MG529

* 3 dilutions of oligos were made :

1:10 1:100 1:1000

1) * ligation rxn. setup as follows :

1 μ l cut plasmid (pEGFP-C1 cut w/ BspEI & BglII)

1 μ l oligo dil'n (1:10 ; 1:100 ; 1:1000)

1 μ l T4 DNA ligase (400000 u/ml)

1 μ l ligase buffer (10x)

6 μ l ddH₂O

10 μ l tot. vol.

* Control : 1 μ l cut plasmid + 9 μ l ddH₂O \Rightarrow 10 μ l
tot. vol.

2) Incubated @ 14-15° overnight

Transformation of DM1 cells w/ pEGFP-C1 ligation rxn.

7-13-96

1) DM1 cells thawed on ice (~100 μ l cells) and ligation rxns. added to cells ; incubate in ice for 30' ; then heat shock cells

2) incubate @ 37° for 1 hr. ; & plated on LB + kan plates which were incubated @ 37° overnight

* 100 μ l were plated on each plate

* very low transformation (1-2 colonies in 1:100 & 1:1000 diln. ; no colonies in 1:10 dil'n) ; will do a new ligation but use HB101 cells instead.

Ligation pEGFP-C1 (cut) w/ oligos

7-15-96

- 1) ligation rxn. setup as before and incubated @ ~~100~~ RT for ~ 7 hrs.
- 2) ligation rxn. were then transformed into HB101 cells; 100 μ l ea. sample plated on kan plates & incubated overnight @ 37°

7-16-96

- 3) started culture of cell colonies in 2 ml LB + kan and incubated 37° overnight
* culture will be made from colonies w/ 1:10 diln of ligation rxn.

7-17-96

- 1) Miniprep of cultured cells containing ligation of pEGFP-C1 and M6528/M6529 oligos;

DNA resuspended in 30 μ l 1XTE

- 2) DNA digested w/ HincII for 2 1/2 hrs. @ 37° in incubator; digestion rxn. was setup as follows:

5 μ l DNA (pEGFP-C1)

0.1 μ l RNase

0.5 μ l HincII

1 μ l Buffer 3

0.05 μ l BSA (mistake!)

3.4 μ l ddH₂O

10 μ l tot. volume

should
be 0.1 μ l

←

30X

3 μ l RNase

15 μ l HincII

30 μ l Buffer 3

1.5 μ l BSA (100x)


100.5 μ l ddH₂O

150 μ l tot. vol.

HinCII digestion of pEGFP-C1 minis

7-17-96

3) DNA run on 2% agarose gel @ 70-100 V :



* used 20 μ l marker ; loaded 12 μ l spl.

* see 147 bp band but no 267 bp or 302 bp bands
since now working w/ pEGFP-C1 parent vector and
not pS65T-C1 ; will react w/ other enzymes

Digestion of pEGFP-C1 (minis)

1) digestion rxns. were setup as follows :

5 μ l DNA

0.1 μ l RNase A

0.1 μ l BSA

0.5 μ l Nhe I

0.5 μ l Xho I

1 μ l Buffer 2

2.8 μ l ddH₂O

10 μ l tot. vol.

30x
3 μ l
3 μ l
15 μ l
15 μ l
30 μ l
84 μ l
150 μ l

2) ~~incubated~~ @ 37° in PCR for 4 hr.

(expected \Rightarrow ~~750~~ bp-insert ; ~~650~~ bp control)


750
..onc

aa

XhoI & NheI digestion of pEGFP-C1 ligation

7-18-96

1) digestion rxn. run on 1.5% agarose gel @ 100V :



* 20 μ l marker ; 12.5 μ l spl. loaded ; overloaded control (contains 5 μ l pEGFP-C1 prep ; next time only use 1 μ l diluted to 5 μ l)

* col. 2 & 4 appear slightly larger than control but cannot say definitively since resolution of 50 bp increase in correct construct is hard to resolve

* will passage 2 & 4 through DM1 and when these colonies are cut w/ Age I & Bsp EI will also cut parent vector w/ same enzymes for a control ; again, correct fragment should be ~50 bp larger than control fragment

Transformation of DM1 w/ Suspected pEGFP-C1 ligation

7-18-96

positives

- 1) transform DM1 cells w/ DNA from colonies 2 and 4;
heat shock @ 42°C for 45 sec.; incubate in LB
@ 37° for 1 hr. and plated on Kan plates (100 μl cells)

pEGFP-C1 ligation & digestion

8-5-96

- * Cultured colonies 2 & 4 in Kan + LB 3ml
- * Miniprep of culture 4 and digestion w/ Age I

* pEGFP-C1 ligation rxn. was digested w/ Age I as follows:

	<u>DNA</u>	<u>Buffer I</u>	<u>Age I</u>	<u>RNAse A</u>	<u>ddH₂O</u>
1) pEGFP-C1 lig	40 μl	10 μl	4 μl	1 μl	45 μl
2) pEGFP-C1 (cont)	1 μl	1 μl	1 μl	—	7 μl

* digested overnight @ RT

8-6-96

* EtOH ppt. Age I digest'n rxn. & digested w/ Bsp EI for 4 hr. @ 37°

* minipreps done on colony 2 followed by Age I digestion as described above

→ 2 control rxns. setup as above

* for Bsp EI digestion, ligation rxn. digested as follows:

	<u>DNA</u>	<u>Buffer I</u>	<u>Bsp I</u>	<u>RNAse A</u>	<u>ddH₂O</u>
1) pEGFP-C1 lig	20 μl	5 μl	4 μl	1 μl	20 μl
2) pEGFP-C1 (control)	10 μl	5 μl	1 μl	—	34 μl

@ 37° for 4 hr.

8-11-96

* started liquid culture of pE6FP-N1 in 3ml LB+amp

8-12-96

* Miniprep 1ml pEGFP-N and digested as follows:

5 μ l DNA (pEGFP-N)	} 30X	3 μ l
0.1 μ l RNase A		3 μ l
1 μ l EcoRI		30 μ l
1 μ l EcoRI Buffer		30 μ l
2.9 μ l ddH ₂ O		87 μ l
10 μ l tot. vol.		

Digestion of Colony 2 w/ BspEI

8-7-96


* digestion rxn. analyzed on 1.5% agarose gel (100V):



DATE: 226 5-77990.60
ID# 565-1441 FILE: (UNTITLED)

* negative ☹

* next time load less
of ligation; can
always run a second in
2 lanes



DATE: 226 5-77990.60
ID# 565-1441 FILE: (UNTITLED)

→ pE6FP
N1

* Same gel allowed to run for
longer time to improve resolution
* though control has very weak
signal they can still be discerned
to be smaller than ligh' rxn.
* 750 band isolated; purified w/
Qiagen column; resuspend in
50 µl TE

Ligation of pFOXSB & pEGFP-N 750 fragment


8-8-96

<u>Label</u>	<u>Contents</u>					
	ddH ₂ O	BspEI cut pFOXSB	pEGFP-N frag 750 bp	ligase	lig. buffer	
	2ul	5ul	1ul	2ul	1ul	1ul
	5ul	2ul	↓	5ul	↓	↓
	5ul	2ul	↓	↓	↓	↓
	10ul	1.5ul	↓	10ul	↓	1.5ul
(control)	C	9ul	—	—	—	—

* rxn. run @ RT for 5 hr.

* transformed into HB101 cells & plated on amp^r

EtOH ppt
 * miniprep & digestion w/ Bsp EI of colony 4 (from yesterday) & will analyze if positive (has 750 band) on 1.5% agarose gel :



* DNA lost during ppt. !

8-11-96

* started liquid culture of pE6FP-N1 in 3ml LB+amp

8-12-96

* miniprep 1ml pE6FP-N and digested as follows:

5ul DNA (pE6FP-N)	} 30X	3ul
0.1ul RNase A		35ul
1ul EcoRI		30ul
1ul EcoRI Buffer		87ul
2.9ul ddH ₂ O		
10ul tot. vol.		


* incubate @ 37° for 4 hrs.

* DNA used was resuspended in 30ul 1X TE

* for control, used 1ul pFoxSB DNA (2 controls made)

* for culture used to make DNA for this digestion, note that tubes 1, 2 & 17 were clean (i.e. apparently, no bacteria) and so should be negatives on gel

* digestion run on 0.5% agarose gel:



* 7 definite positives

(i.e. has 850 & 175 band)

5, 7, 9-11, 20 & 22;

possible positives include

4 & 13

→ will be

pE6FP-N1

8-13-96

* pEGFP-N1 from miniprep digested w/ Sac I :

27 μ l DNA (#10 pEGFP-N1)

1 μ l BSA

1 μ l RNase A

4 μ l Sac I

5 μ l Buffer I

12 μ l ddH₂O

50 μ l

* digested @ 37° for 4 hr.

* reculture #5 & #10 into ~~new~~ fresh LB+amp (3 ml)

* pEGFPN1-SacI digest CIP 2x @ 37° for 30' w/ 1 μ l CIP

in following rxn:

50 μ l SacI digest

30 μ l 10x CIP buffer

1 μ l CIP

219 μ l ddH₂O

300 μ l tot. vol.

* CIP rxn. Stop by adding: 39 μ l ddH₂O

and heating @ 68° for 15'

40 μ l 10x STE

20 μ l 10% SDS

400 μ l tot. vol. w/ SacI digest

* extract w/ phenol/chloroform 2x; EtOH ppt & resuspend

in 15 μ l TE (1x)

* run on 1% agarose gel & ~~purified for ligation rxn:~~

to assay concentration (1 μ l DNA assayed) w/

6 μ l marker (1x)


8-13-96

* Assay for concentration of pEGFP-N1 cut SacI & CIP
 on 1% agarose gel ; 6 μ l marker (1x) :

* ~20 ng/ μ l pEGFP-N1 to be used
 for ligation

* Ligation reactions setup as follows:

Name	PEGFP-N1 cut SacI/CIP	insert control	Ligase Buffer (10x)	Ligase 400000u	ddH ₂ O
2 μ l	3 μ l	2 μ l	1 μ l	1 μ l	3 μ l
6 μ l	↓	6 μ l	1.5 μ l	↓	3.5 μ l
control	↓	—	—	—	7 μ l



* ligated @ 15° overnight

8-14-96 * transformed into 100 μ l HB101 cells & plated 100 μ l onto LB+amp
 @ 37° overnight

* Qiagen column (20-tip) on colonies 5 & 10 (pEGFP-N1)
 and resuspended in 20 μ l 1x TE

* Digested pF0CAT2 w/ SacI to isolate 175 bp β -globin
 intron:

50 μ l PFC2

4 μ l SacI

1 μ l BSA

7.5 μ l Buffer I


12.5 μ l ddH₂O

75 μ l tot. vol.

@ 37° for 4 hr.

8-14-96

* pEC2 SacI digest to isolate 175 bp intron
run on 1.5% gel :



* isolated & purified 175 bp β -globin
intron w/ Qiagen ; resuspended in
50 μ l 1x TE

303 34175055 04/09/96 11:55
MEZ 18.12 ID# 669-01546 FILED (UNTITLED)

8-15-96

* plates of pEGFP-N1 ligation w/ 175 bp intron show hi
transformation efficiency ; will culture 12 colonies
each (w/ total of 24 cultures) from 2 μ l & 6 μ l
plate into 3 ml LB + amp ; 37 $^{\circ}$ overnight

* cut #5 & #10 pEGFP-N1 w/ SacI and kept in
freezer ; digest 10 μ l ea. DNA w/ 4 μ l enzyme in
50 μ l rxn.


8-16-96

* minis of suspected pEGFP-N2 ; minis on 1 ml culture ;
after add'n of 900 μ l 100% EtOH put in EtOH bath
for 15' and centrifuged for 10' w/ subsequent
aspiration ; resuspended in 30 μ l 1X TE

* susp. pEGFP-N2 digest w/ EcoRI

5 μ l DNA (susp. pEGFP-N2 ; pEGFP-N1 control ⁺ used 1 μ l)		
0.1 μ l RNase A	} 35x	
1 μ l EcoRI		
1 μ l EcoRI Buffer		
2.9 μ l ddH ₂ O		
10 μ l tot. vol.		
		3.5 μ l RNase A
		35 μ l EcoRI
		35 μ l EcoRI Buffer
		101.5 μ l ddH ₂ O
		175 μ l

* digest run on 1.5% gel :




EXPOSURE: 6/30 SEC. E-BOX 0 IMAGE 255 BYTIME 1.05
DATE: 08-23-1996 TIME: 18:44 ID# 569-0153 FILE: (UNTITLED)

* recultured #4, 6, 19, 22,
23 in

3ml LBamp

* QiaGen column on
#4, 6



#4
~180 ng
#6
~60 ng

SEC. E-BOX 0 IMAGE 255 BYTIME 1.05
DATE: 08-23-1996 TIME: 18:44 ID# 569-0153 FILE: (UNTITLED)

8-19-96

* pE6FP-N2 digested again w/ EcoRI & set up
dbl digestion w/ EcoRI & BspMI (will repppt.)

5 μ l DNA

1 μ l EcoRI (BspMI)

1 μ l EcoRI Buffer (Buffer 2)

3 μ l ddH₂O

10 μ l tot. vol.

* digested pFC2 - 482 rglcZ w/ HindIII & BglII

30 μ l DNA

5 μ l Buffer B (Boehringer Buffer)

2.5 μ l HindIII

2.5 μ l BglII

10 μ l ddH₂O

50 μ l tot. vol.

all @ 37° 4 hr.

* second set dbl digestion on # 4, 6, 8, 12, 14, 15, 21-23

15 μ l DNA

1 μ l RNase A

2 μ l BspMI

10 μ l Buffer 2

72 μ l ddH₂O


100 μ l tot. vol.

incubate 37° overnight

10x

* digested pFC2 - 482 rglcZ run on 1% agarose
and purified ~730 bp fragment:

* resuspended in 50 μ l TE



EXPOSURE: 1.0 SEC. PLATE: 255
DATE: 08-19-1996 TIME: 15:19 TDO 669-016

8-20-96

* Digestion w/ Bsp MI & Eco RI to check for intron orientation:

* yesterday's Bsp MI rxn. phenol / chloroform extracted (1x)
ppt. w/ cold EtOH; resuspended in 10 μ l 1x TE;
digested w/ Eco RI as follows:


10 μ l DNA
1 μ l Eco RI
2 μ l buffer
0.1 μ l RNase A
6.9 μ l ddH₂O
20 μ l tot. vol

} 15x

15 μ l Eco RI
20 μ l buffer
1.5 μ l RNase A
103.5 μ l ddH₂O

* digested 4 hr 37°

* run on 1.5% agarose gel; culture #4 & 23




* can see 905 bp band on #4 (tip) and can be discerned for #8, 22, 23; bad digestion w/ Bsp MI since can't see 250 bp band for negatives # 12, 14, 15, 21

EXPOSURE: 8.0 SEC. BLACK: 0.75 WHITE: 255.0 DENSITY: 0.75 DATE: 08-20-1996 TIME: 12:54 FILE: 569-0.534 FILE: (UNITED)

8-21-96

* digested pEGFP-N2 w/ Hind III & Bgl II so com
 ligate pFC2 -482 rglc Z :
 10 μ l DNA (~10 μ g)
 2.5 μ l Hind III
 2.5 μ l Bgl II
 5 μ l buffer B
 30 μ l ddH₂O
 50 μ l tot. vol. @ 37° for 4 hr.



* run on 0.8% agarose gel ;
 Qiagen gel extract vector & resuspend
 in 50 μ l

* setup ligations w/ -482 rglc Z 73C
 fragments

cutEGFP	insert	ligase	lig. buffer	ddH ₂ O
1 μ l	2 μ l	1 μ l	1 μ l	5 μ l
↓	5 μ l	↓	1.5 μ l	6.5 μ l

* control = 1 μ l cutEGFP + 9 μ l ddH₂O

* ligations incubate @ 15° overnight

8-22-96

* plated ligation rxns onto amp + LB plates (100 μ l ea.)

8-23-96


* started 24 cultures in ~3 ml LB + amp

* picked from "2 μ l insert" plate

* did not observe in transformation efficiency

.10/1. mmis on cmv. 40 kip. rgl. w GFP. V W RE (Janet)







Sample 08
DyeTerminator(AnyPrimer)
Lane 4

Signal: G:62 A:71 T:37 C:22


Points 60 6184 Base 1: 604
Instrument #907443
FD 2318

Wed, Aug 21, 1996 2:59 PM Page 1 of 2
X: 0 to 6947 Y: 0 to 1200
Spacing: 10.51





Spacing: 10.91





Model 373A

Version 1.2.0

GTACCCCAAGACATNAAGCNCACAAATGCTTCATTCT

Sample 19
Dye Terminator(AnyPrimer)
Lane 19


Signal: G:352 A:485 T:305 C:158


Points 69
Instrument #907443
Control pgEM

Base 1: 691

Wed, Aug 21, 1996 2:59 PM
X: 0 to 6503 Y: 0 to 1200
Spacing: 11:03

Page 2







Model 373A

Version 1.2.0


Sample 18.1
Dye Terminator (Any Primer)
Lane 18

Signal: G:358 A:486 T:319 C:157

5' TCCAGCCGC TAC CCC ACC AGA TNAANCAAG ACAATTCTTGGCANTCT
480 586 610

Points 606 6184 Base 1: 690
Instrument #907443
FD 2326

Wed, Aug 21, 1996 2:59 PM Page 2 of 2
X: 0 to 6631 Y: 0 to 1200
Spacing: 10.82



8-15-90

* Maize's frag. → used 50ng vector
~~used 1.0 insert~~ 3.0 insert

* 482 rgs promoter (rat glucagon)

* pEGFP-N2 = do double digest only w/ pos.

carrier DNA: digest ; when ppt. w/ EtOH add 1.0 μl
tRNA then continue

* plasmid prep → dbe CsCl (2x) of pEGFP-N2

* also reculture & do 20-40 tip of pEGFP-N2

* always leave 1 μl of DNA from prep.

* transform 1 μl DNA for a plate

~~for~~

* for BamHI / HindIII digest use Buffer B ;
can phenol / chloroform purify or gel purify

* intron - inc. exp. level & thus sens. of system


*


8-24-96

* minus pEGFP-N2 -482 rglc Z- element w/ EcoRI
overnight

pick minus for pEGFPN2-482+rglc 2 element

9-25-9





EXPOSURE: 12/30 5Z. BLACK: 0 WHITE: 255 GAMMA: 1.05
DATE: 08-25-1998 TIME: 09:56 ID: 569-01698 FILE: (UNTITLED)