

7-1-96

MIDI Prep of Plasmid DNA from Transformed DM1

- 1) culture started from colony w/ 100 μ l cells (see previous) plated @ 3 μ ml \rightarrow then 30 ml
- 2) MIDI Prep of plasmid DNA (see Qiagen MIDI Prep protocol)
- 3) DNA resuspended in \sim 30 μ l TE

Notes: * culturing of DM1 was very slow; about 8 hr. to get from clear media to turbid
* centrifugation to pellet DNA during MIDI prep was not favorable; DNA not readily visible so DNA ~~may~~ (or sufficient quantities of DNA) may not have been isolated

7-1-96 → 7-2-96

Digestion of pEGFP-N w/ Age-I & Bsp EI and pFOXSB w/ Bsp EI

total rxn-
volumes100 μ l7-2-96

- 1) pEGFP-N (³⁰~~25~~ μ g from previous experiment) was digested w/ Age I @ room temp. overnight (4 μ l enzyme)
- 2) DNA was reprecipitated w/ cold EtOH, resuspended and cut w/ Bsp EI (4 μ l) for ~4 hr. @ 37°C
- 3) Concurrently, pFOXSB (5 μ g) was cut w/ Bsp EI (4 μ l) and incubated @ 37° for ~4 hr.
- * pEGFP-N (cut w/ Bsp EI) and pFOXSB was digested in 50 μ l total volume rxn.
- 4) Digested pFOXSB was then mixed w/ CIP for ²⁵~~20~~ minutes

- 5) Samples run on 1.5% agarose gel (100 V) w/ λ marker



* used only
1/2 amt. of
EtBr in
agarose

pFOXSB

pEGFP-N

 λ

cut w/:

Bsp EI

Bsp EI

+ Age I

7-2-96

6) pFOXSB (BspEI cut) band and pEGFP-N 750 bp band was cut out from gel and purified using Qiagen Gel Extraction Kit

* pFOXSB band = 215 mg

* pEGFP-N 750 bp band = 206 mg

7) DNA eluted from spin column was resuspended in 50 μ l 1x TE

Ligation of pFOXSB and pEGFP-N 750 bp fragment (EGFP + NLS)

1) The following ligation reactions were setup:

<u>Tube #</u>	<u>contents</u>
1	1 μ l cut pFOXSB + 5 μ l pEGFP-N fragment + 1 μ l ligase + 1 μ l ligase buffer + 2 μ l ddH ₂ O (ligase added last) tot. vol = 10 μ l
2	1 μ l cut pFOXSB + 10 μ l pEGFP-N fragment + 1 μ l ligase + 1.5 μ l buffer + 1.5 μ l ddH ₂ O tot. vol = 15 μ l
(control) 3	1 μ l cut pFOXSB + 9 μ l ddH ₂ O (for background) tot. vol = 10 μ l

2) The above contents were mixed (except for ligase) in Eppi & centrifuged briefly; then ligase was added

Ligation (cont'd)

7-2-96

3) ligation rxn. incubated overnight
@ 14-15°C

7-3-96

Plasmid MIDI Prep

* pellet
bacteria
from culture
by spinning
@ 40° for 10'
speed 4

1) MIDI Prep performed on DM1 transformed cells from a second culture (prepared by Maibe) following culture procedures described in Qiagen plasmid MIDI protocol. (see 7-1-96)
* DM1 culture yielded more cells this time (~3X)

2) for MIDI prep., P1 buffer and P2 buffers were added before resuspension ~~and mixing~~ of cells.

3) After ~~15' incubation in ice~~ the addition of P3 buffer and subsequent incubation on ice for 15', the mixture was aliquoted into Eppi's w/ 1 ml in each tube for a total of 12 tubes.

4) The tubes were then centrifuged @ hi speed for 30' ~~in microcentrifuge~~ in microcentrifuge (@ RT (room temp.))

5) Qiagen tip-100 was equilibrated w/ 4 ml QBT buffer

6) Supernatant from each tube was applied to column

7) Column washed 2x10 ml QC buffer

8) DNA eluted from column w/ QF buffer (5ml) and collected in Eppi's (aliquoted to 500 µl per tube)

9) DNA was ppt. by adding 350 µl isopropanol to each tube & centrifuging @ hi speed (RT) for 30'


Plasmid MIDI Prep (cont'd)

7-3-96

10) Ethanol (cold & 70%) was used to wash DNA (500 μ l / tube) and redissolved in 75 μ l TE 1X
 * 75 μ l TE was added to tube 1 and this same TE was used to resuspend ppt. DNA in the rest of the 9 tubes; ultimately, all DNA (originally separated in 10 tubes) is resuspended in 1 tube labelled "pEGFP-N1 (MIDI) 7-3-96"

run @
100V

11) an aliquot of this prep was then run on an agarose (1%) gel:



pGK-neo (Midi)
pEGFP-N1

* 10 μ l λ marker

* pEGFP-N1 & pGK-neo = 1 μ l plasmid

7-3-96

Transformation of Competent cells w/ pFOXEGFP-N1 ligation

*bacteria
used is
HB101

- 1) cells thawed on ice (100ul cells/tube), then add all of ligation rxns. to tubes of cells
- 2) a fourth transformation rxn. was performed using BS ΔNot to transform cells (1ul BS ΔNot) for Maïke
- 3) incubate on ice for 30'; warm 0.9 ml LB @ 42°
- 4) heat shock @ 42° for 45 sec and ice for 2'; then add ligation rxns. to ~~cells~~ ^{LB} (all of the rxn. tube)
- 5) incubate @ 37° for 1 hr.; briefly spin down cells (5 sec.)
- 6) plate on amp plates as follows:
 - a) plate 100ul of 5ul insert, 10ul insert and control
 - b) streak plate w/ loop 5ul BS ΔNot

* 1ul BS ΔNot → 10ul total vol. w/ ddH₂O; then added to 100ul competent cells
- 7) incubate @ 37° ~~for 1 hr.~~ ^{remove 900ul & resuspend cells in remaining 100ul media} overnight
- 8) picked 24 colonies from plate w/ 10ul insert and grown in culture tubes w/ 3ml LB overnight in shaker

7-5-96

Miniprep of pFOXEGFP-N1 and Endonuclease Digestion

- 1) put 1 ml of liquid culture (from total of 24 cultures) in Eppi's and spin 50 sec.
- 2) pour off supernatant by ~~flicking~~ inverting tube once; this leaves ~50 μ l in tube; resuspend pellet in remaining media by vortexing
- 3) add 300 μ l TENS (cell lysing) ^{mix} and 150 μ l 3M NaOAc pH 5.2 ^{mix}; centrifuged @ hi speed for 3'
- 4) supernatant transferred to new Eppi's and added 900 μ l 100% cold EtOH; mix by inverting tubes; then, centrifuge @ hi speed for 3'
- 5) aspirate off supernatant and wash w/ cold 70% EtOH; ~~gently invert tubes~~ and centrifuge @ hi speed for 5' (do not mix or invert tubes!)
- a) aspirate off supernatant and resuspended in 30 μ l 1x TE

Digestion

w/ EcoRI

- 1) digestion of pFOXEGFP-N1 w/ EcoRI were setup as follows:


master mix prepared on ice	{	5 μ l DNA (pFOXEGFP-N1; and pFOXSB as control)
		0.1 μ l RNase (10 mg/ml)
		0.5 μ l EcoRI
		1 μ l Buffer I
		3.4 μ l ddH ₂ O
		10 μ l total volume

4

Digestion of pFOXEGFP-N1 (cont'd)

7-5-96

- 2) digestion @ 37° was performed for 2 hrs.
- 3) add 2 μ l loading dye (6X) to each tube and load samples on 1.5% agarose gel and run @ 90 V :



- * 2 and 13 = no DNA
- * something wrong w/ gel?
- * DNA does not appear digested well ; digest > 2hr?
- * not well digested b/c used Buffer 1 V
- * control = pFOXSB (1 μ l plasmid)

7-8-96

Reculture of colony 24 for miniprep

- 1) put sterile ~~amp~~ LB (3 ml) into 2 culture tubes each (LB + amp)
- 2) flame an inoculating loop, cool and dip into tube w/ colony 24 and dip into fresh tube w/ LB to reculture; repeat a second time and incubate tubes @ 37° overnight until turbid;
additional tubes (for Maize):

② BS Δ Not

⑤ pGK-neo

* for these, inoculated a colony from each plate and cultured overnight in 3 ml LB + amp;

Miniprep (Qiagen) of colony 24 & BS Δ Not &

7-9-96

pGK-neo

* bacteria
used is
HB101

- ① pellet bacteria by spinning (speed 4) @ 4° C for 10'; discard supernatant
- ② add 300 μl Buffer P1 and resuspend pellet; transfer to Eppi's
- ③ add 300 μl Buffer P2 gently mix and incubate @ RT for 5'
- ④ add chilled Buffer P3 (300 μl); invert 5-6 x and incubate on ice for 10'
- ⑤ centrifuged @ hi speed for 15'; supernatant moved to 20-tip column (equilibrated)
- ⑥ wash 4 x 1 ml Buffer QC, then DNA eluted w/ 800 μl Buffer QF

Miniprep (cont'd)

7-9-96

⑦ DNA ppt. w/ 560 μ l isopropanol (0.7 vol)
and centrifuged @ hi speed for 30'

⑧ DNA washed w/ 1 ml 70% cold EtOH;
centrifuge 15 sec. and redissolved in 20 μ l 1x TE

Digestion of Colony 24

1) Digestion set up as follows: (2 rxns. made)

1 μ l DNA (pFOXEGFP-N1)

1 μ l EcoRI Buffer

1 μ l EcoRI

7 μ l ddH₂O


10 μ l tot. vol.

Digested for 1 hr. @ 37°C

2) loaded onto 1.5% agarose gel:

* 20 μ l marker

* 12 μ l digestion sp



* put in
EtBr
bath
1 μ g/ml
for 10';
wash for
5'

* where's the 850 bp band?

* funky gel

* ~ 1000 bp band seen


* will redigest


* 10 weight band


observed on

minitor. print

AATCTGTTGAA









Model 373A

Version 1.2.0

Sample 12

Dye Terminator (Any Primer)

Lane 12

Signal: G:409 A:331 T:257 C:164

Points 62

Instrument: 807443

MS 2170

Base 1: 620


Wed, Jul 10, 1996 3:00 PM

X: 0 to 6867 Y: 0 to 1200

Spacing: 10.48

Page 2 of 2

5' JNNNAC GNA 486 T GCA AAC AAN C T NNGGGA AA AA A C T CCA 830 T AACTCCA 840 C AANCC NC



Sac I Digest pFOXEGFP-N1

7-9-96

1) the following digest was setup:

10 μ l DNA (pFOXEGFP-N1)

1 μ l BSA (final conc. 2X)

4 μ l Sac I

5 μ l Buffer I

30 μ l ddH₂O


50 μ l total volume

* a second digestion w/ Eco RI was setup as previously described (1 μ l Eco RI & 1 μ l DNA); both digestions performed in PCR machine @ 37° for 4 hr.

7-10-96

Eco RI digest of pFOXEGFP-N1 was run on 1.5% agarose gel:

* colony
24 #1



* 20 μ l λ marker & 12 μ l rxn. Spl

* bad digestion: can see a lot of undigested plasmid; also observe the partial digest band ~ > 1000 bp

→ this may have resulted from incomplete digest to make 850 + 175 \Rightarrow 1025 b band

* 850 bp (GFP & NLS) band observed

CIP of pFOXEGFP-N1

7-10-96

- ① CIP rxn. of Sac I digested vector was setup as follows: (added directly to digestion rxn.)

30 μ l 10X CIP Buffer

219 μ l ddH₂O

1 μ l CIP

+ 50 μ l Sac I digest

300 μ l total vol. ; incubate @ 37° for 30'

- ② add 1 μ l additional CIP, then incubate @ 37° for 30'

- ③ stop rxn. by adding :

(39) 40 μ l ddH₂O

40 μ l 10X STE

20 μ l 10% SDS

400 μ l total volume ; heat @ 68° for 15'

- ④ Extract 2X w/ phenol / chloroform in 1:1 (v/v) ratio

* add phenol / chloroform; vortex; centrifuge @ hi speed for 5'

- ⑤ EtOH precipitate DNA w/ 1 ml cold EtOH ; incubate in dry ice / EtOH bath for 15' ; centrifuge 15' ; wash w/ 1 ml 70% EtOH ; resuspend in 50 μ l 1X TE

15

Ligation of pFOXEGFP-N1 w/ human β -globin intron

7-10-96

1) ligation rxns. were setup as follows:

- a) 1 μ l plasmid (pFOXEGFP-N1)
2 μ l insert (β -globin intron)
1 μ l ligase ^{175 bp} (400000 u/ml)
1 μ l ligase buffer (10x)
5 μ l ddH₂O
10 μ l tot. vol.

- b) 1 μ l plasmid
6 μ l insert
1 μ l ligase
1.5 μ l lig. buffer
5.5 μ l ddH₂O
15 μ l tot. vol.

- c) control: 1 μ l plasmid + 9 μ l ddH₂O = 10 μ l tot. vol.

* 175 bp β -globin human intron was not quantified
2) ligation rxns. were incubated @ 16° overnight

Transformation of HB101 cells w/ pFOXEGFP-N1

7-11-96

Ligation

1) cells thawed on ice and added directly to ligation rxns.; rest of transformation performed as previously described; plate 100 μ l rxn. to each plate (LB + amp)

2) plate on LB+amp (100 μ l) of

- a) control b) 2 μ l insert c) 6 μ l insert

* incubate overnight (start 12:00 \Rightarrow 4pm)

7-11-96

7-12-96

Miniprep of pFOXSB

7-11-96

- 1) as described in Sbera's protocol, except put on ice for ~~5~~ 5' after adding NaOAc (3M pH 5.2)
- 2) after adding 100% EtOH, put in methanol/CO₂ bath for 15'; centrifuge @ hi speed for 3'
- 3) DNA resuspended in 30 μ l 1x TE

Neo I Digestion of primary plasmids used to construct GFP plasmid

* sequence data found on next page

* sequence data of pFOXEGFP-N1 showed that the wrong plasmid may have been used to construct GFP plasmid; cutting w/ NcoI should determine whether the right plasmid (pEGFP-C1) and not the wrong old plasmid (pSG5T-C1) was used to construct GFP plasmid. The wrong plasmid pSG5T-C1 will yield an extra low weight fragment (167 bp).

1) Digestion rxn. setup as follows:

Mix made

5 μ l DNA
0.1 μ l RNase
1 μ l NcoI
1 μ l Buffer 4
2.9 μ l ddH₂O
10 μ l tot. vol.

* for commercial plasmids pEGFP-C1 and pSG5T-C1 only 1 μ l DNA was used and added to 4 μ l ddH₂O; and this was added to 5 μ l of mix w/ a tot. rxn. vol. of 10 μ l


* also digested minis (pEGFP-N) # 37-48 (positives) @ 37° for 1 1/2 hrs.; then stored in freezer

* digestion of pEGFP-C1 (prep); pEGFP-N (MIDI prep) and pEGFP-N (unmethylated) was performed

~~Not I~~ NcoI Digestion (cont'd)

7-12-96

i) digestion rxns. were run on 2% agarose gel : (120V)



same order:

- | | | |
|--------------------------------------|--------------------------|-------------------|
| 1) Marker | 2) pEGFP-C1 (commercial) | 3) pS65T-C1 (com) |
| 4) pEGFP-C1 (prep) | 5) pEGFP-N (methylated) | |
| 6) pEGFP-N (MIDI prep, unmethylated) | | |

* used 20ul marker

* ~ 12 μ l sp loaded

→ presence of extra 167 bp fragment indicates minis made from old parent vector pSG5T-C1

→ start over w/ pEGFP-C1 prep

7-12-96

Annealing of Oligos

M6528 &
M6529 @
conc. of
1 $\mu\text{g}/\mu\text{l}$

1) added 1 μg (1 μl) of oligo M6528 and M6529 into rxn. w/ 2.5 μl Buffer H (Boehringer-Mannheim) to a total rxn. vol. of 25 μl

2) boiled in water bath for ~1 minute then kept in water bath as it eq. temp. equilibrated to RT

* no kinase

3) put @ 4°C for 2 hr. before using
Digestion of pEGFP-C1 (prep)

1) digestion rxn. setup as follows:

10 μl pEGFP-C1 (prep)

2.5 μl Bsp E1

2.5 μl Bgl II

5 μl Buffer 3

30 μl ddH₂O

50 μl total volume

2) incubate @ 37° for 4 hr.

3) phenol / chloroform and ethanol ppt. the DNA ;
resuspend in 30 μl 1xTE

Ligation of pEGFP-C1 (cut) w/ oligos

7-12-96

MG528 / MG529

* 3 dilutions of oligos were made :

1:10 1:100 1:1000

1) * ligation rxn. setup as follows :

1 μ l cut plasmid (pEGFP-C1 cut w/ BspEI & BglII)

1 μ l oligo dil'n (1:10 ; 1:100 ; 1:1000)

1 μ l TH DNA ligase (400000 u/ml)

1 μ l ligase buffer (10x)

6 μ l ddH₂O

10 μ l tot. vol.

* Control : 1 μ l cut plasmid + 9 μ l ddH₂O \Rightarrow 10 μ l
tot. vol.

2) Incubated @ 14-15° overnight

Transformation of DM1 cells w/ pEGFP-C1 ligation rxn.

7-13-96

1) DM1 cells thawed on ice (~100 μ l cells) and ligation rxns. added to cells ; incubate in ice for 30' ; then heat shock cells

2) incubate @ 37° for 1 hr. ; & plated on LB + kan plates which were incubated @ 37° overnight

* 100 μ l were plated on each plate

* very low transformation (1-2 colonies in 1:100 & 1:1000 diln. ; no colonies in 1:10 dil'n) ; will do a new ligation but use HB101 cells instead.

Ligation pEGFP-C1 (cut) w/ oligos

7-15-96

- 1) ligation rxn. setup as before and incubated @ ~~100~~ RT for ~ 7 hrs.
- 2) ligation rxn. were then transformed into HB101 cells; 100 μ l ea. sample plated on kan plates & incubated overnight @ 37°

7-16-96

- 3) started culture of cell colonies in 2 ml LB + kan and incubated 37° overnight
* culture will be made from colonies w/ 1:10 diln of ligation rxn.

7-17-96

- 1) Miniprep of cultured cells containing ligation of pEGFP-C1 and M6528/M6529 oligos;

DNA resuspended in 30 μ l 1XTE

- 2) DNA digested w/ HincII for 2 1/2 hrs. @ 37° in incubator; digestion rxn. was setup as follows:

5 μ l DNA (pEGFP-C1)

0.1 μ l RNase

0.5 μ l HincII

1 μ l Buffer 3

0.05 μ l BSA (mistake!)

3.4 μ l ddH₂O

10 μ l tot. volume

} 30X

3 μ l RNase

15 μ l HincII

30 μ l Buffer 3

1.5 μ l BSA (100x)

100.5 μ l ddH₂O


150 μ l tot. vol.

should be 0.1 μ l

HinCII digestion of pEGFP-C1 minis

7-17-96

3) DNA run on 2% agarose gel @ 70-100 V :



* used 20 μ l marker ; loaded 12 μ l spl.

* see 147 bp band but no 267 bp or 302 bp bands
since now working w/ pEGFP-C1 parent vector and
not pSG5T-C1 ; will react w/ other enzymes

Digestion of pEGFP-C1 (minis)

1) digestion rxns. were setup as follows :

5 μ l DNA

0.1 μ l RNase A

0.1 μ l BSA

0.5 μ l Nhe I

0.5 μ l Xho I

1 μ l Buffer 2

2.8 μ l ddH₂O

10 μ l tot. vol.

3 μ l
3 μ l
15 μ l
15 μ l
30 μ l
84 μ l
150 μ l

30x

2) ~~incubated~~ @ 37° in PCR for 4 hr.

(expected \Rightarrow ~~750~~ bp-insert ; ~~650~~ bp control)

750


..onc

aa

XhoI & NheI digestion of pEGFP-C1 ligation

7-18-96

1) digestion rxn. run on 1.5% agarose gel @ 100V :



* 20 μ l marker ; 12.5 μ l spl. loaded ; overloaded control (contains 5 μ l pEGFP-C1 prep ; next time only use 1 μ l diluted to 5 μ l)

* col. 2 & 4 appear slightly larger than control but cannot say definitively since resolution of 50 bp increase in correct construct is hard to resolve

* will passage 2 & 4 through DM1 and when these colonies are cut w/ Age I & Bsp EI will also cut parent vector w/ same enzymes for a control ; again, correct fragment should be ~50 bp larger than control fragment

Transformation of DM1 w/ Suspected pEGFP-C1 ligation

7-18-96

positives

- 1) transform DM1 cells w/ DNA from colonies 2 and 4;
heat shock @ 42°C for 45 sec.; incubate in LB
@ 37° for 1 hr. and plated on Kan plates (100 μl cells)

pEGFP-C1 ligation & digestion

8-5-96

- * Cultured colonies 2 & 4 in Kan + LB 3ml
- * Miniprep of culture 4 and digestion w/ Age I

* pEGFP-C1 ligation rxn. was digested w/ Age I as follows:

	<u>DNA</u>	<u>Buffer I</u>	<u>Age I</u>	<u>RNase A</u>	<u>ddH₂O</u>
1) pEGFP-C1 lig	40 μl	10 μl	4 μl	1 μl	45 μl
2) pEGFP-C1 (cont)	1 μl	1 μl	1 μl	—	7 μl

* digested overnight @ RT

8-6-96

* EtOH ppt. Age I digest'n rxn. & digested w/ Bsp EI for 4 hr. @ 37°

* minipreps done on colony 2 followed by Age I digestion as described above

→ 2 control rxns. setup as above

* for Bsp EI digestion, ligation rxn. digested as follows:

	<u>DNA</u>	<u>Buffer I</u>	<u>Bsp I</u>	<u>RNase A</u>	<u>ddH₂O</u>
1) pEGFP-C1 lig	20 μl	5 μl	4 μl	1 μl	20 μl
2) pEGFP-C1 (control)	10 μl	5 μl	1 μl	—	34 μl

@ 37° for 4 hr.

8-11-96

* started liquid culture of pE6FP-N1 in 3ml LBamp

8-12-96


* Miniprep 1ml pEGFP-N and digested as follows:

5 μ l DNA (pEGFP-N)	} 30X	3 μ l
0.1 μ l RNase A		3 μ l
1 μ l EcoRI		30 μ l
1 μ l EcoRI Buffer		30 μ l
2.9 μ l ddH ₂ O		87 μ l
10 μ l tot. vol.		

Digestion of Colony 2 w/ BspEI

8-7-96


* digestion rxn. analyzed on 1.5% agarose gel (100V):



DATE: 226 5-7700 0.60
ID# 569-1441 FILE: (UNTITLED)

* negative ☹

* next time load less
of ligation; can
always run a second in
2 lanes



DATE: 226 5-7700 0.60
ID# 569-1441 FILE: (UNTITLED)

* Same gel allowed to run for
longer time to improve resolution
* though control has very weak
signal they can still be discerned
to be smaller than ligh' rxn.
* 750 band isolated; purified w/
Qiagen column; resuspend in
50 µl TE

→ DE6FF
N1

Ligation of pFOXSB & pEGFP-N 750 fragment


8-8-96

<u>Label</u>	<u>Contents</u>					
	<u>ddH₂O</u>	<u>BspEI cut pFOXSB</u>	<u>pEGFP-N frag 750 bp</u>	<u>ligase</u>	<u>lig. buffer</u>	
	2ul	5ul	1ul	2ul	1ul	1ul
	5ul	2ul	↓	5ul	↓	↓
	5ul	2ul	↓	↓	↓	↓
	10ul	1.5ul	↓	10ul	↓	1.5ul
(control)	C	9ul	—	—	—	—

* rxn. run @ RT for 5 hr.

* transformed into HB101 cells & plated on amp^r

* EtOH ppt & digestion w/ Bsp EI of colony 4 (from yesterday) & will analyze if positive (has 750 band) on 1.5% agarose gel :



* DNA lost during ppt. !

8-11-96

* started liquid culture of pE6FP-N1 in 3ml LB+amp

8-12-96

* miniprep 1ml pE6FP-N and digested as follows:

5 μ l DNA (pE6FP-N)	} 30X	3 μ l
0.1 μ l RNase A		30 μ l
1 μ l EcoRI		30 μ l
1 μ l EcoRI Buffer		87 μ l
2.9 μ l ddH ₂ O		
10 μ l tot. vol.		


* incubate @ 37° for 4 hrs.

* DNA used was resuspended in 30 μ l 1X TE

* for control, used 1 μ l pFoxSB DNA (2 controls made)

* for culture used to make DNA for this digestion, note that tubes 1, 2 & 17 were clean (i.e. apparently, no bacteria) and so should be negatives on gel

* digestion run on 0.5% agarose gel:



* 7 definite positives

(i.e. has 850 & 175 bands)

5, 7, 9-11, 20 & 22;

possible positives include

4 & 13

→ will be

pE6FP-N1

8-13-96

* pEGFP-N1 from miniprep digested w/ Sac I :

27 μ l DNA (#10 pEGFP-N1)

1 μ l BSA

1 μ l RNase A

4 μ l Sac I

5 μ l Buffer I

12 μ l ddH₂O

50 μ l

* digested @ 37° for 4 hr.

* reculture #5 & #10 into ~~new~~ fresh LB+amp (3 ml)

* pEGFPN1-SacI digest CIP rx @ 37° for 30' w/ 1 μ l CIP

in following rxn:

50 μ l SacI digest

30 μ l 10x CIP buffer

1 μ l CIP

219 μ l ddH₂O

300 μ l tot. vol.

* CIP rxn. Stop by adding: 39 μ l ddH₂O

and heating @ 68° for 15'

40 μ l 10x STE

20 μ l 10% SDS

400 μ l tot. vol. w/ SacI digest

* extract w/ phenol/chloroform rx ; EtOH ppt & resuspend

in 15 μ l TE (1x)

* run on 1% agarose gel & ~~purified for ligation rxn~~

to assay concentration (1 μ l DNA assayed) w/

6 μ l marker (1x)


8-13-96

* Assay for concentration of pEGFP-N1 cut SacI & CIP
 on 1% agarose gel ; 6 μ l marker (1x) :

* ~20 ng/ μ l pEGFP-N1 to be used
 for ligation

* Ligation reactions setup as follows:

Name	PEGFP-N1 cut SacI / CIP	insert control	Ligase Buffer (10x)	Ligase 400000u	ddH ₂ O
2 μ l	3 μ l	2 μ l	1 μ l	1 μ l	3 μ l
6 μ l		6 μ l	1.5 μ l	\downarrow	3.5 μ l
control	\downarrow	—	—	—	7 μ l



* ligated @ 15° overnight

8-14-96 * transformed into 100 μ l HB101 cells & plated 10 μ l onto LB+amp
 @ 37° overnight

* Qiagen column (20-tip) on colonies 5 & 10 (pEGFP-N1)
 and resuspended in 20 μ l 1x TE

* Digested pFosCAT2 w/ SacI to isolate 175 bp β -globin
 intron :

50 μ l PFC2

4 μ l SacI

1 μ l BSA

7.5 μ l Buffer I


12.5 μ l ddH₂O

75 μ l tot. vol.

@ 37° for 4 hr.

8-14-96

* pEC2 Sac I digest to isolate 175 bp intron
run on 1.5% gel :



* isolated & purified 175 bp β -globin
intron w/ Qiagen ; resuspended in
50 μ l 1x TE

303 3475255 04/04/96 11:55
ME: 18 12 ID# 669-01546 FILE: (UNTITLED)

8-15-96

* plates of pEGFP-N1 ligation w/ 175 bp intron show hi
transformation efficiency ; will culture 12 colonies
each (w/ total of 24 cultures) from 2 μ l & 6 μ l
plate into 3 ml LB + amp ; 37 $^{\circ}$ overnight

* cut #5 & #10 pEGFP-N1 w/ Sac I and kept in
freezer ; digest 10 μ l ea. DNA w/ 4 μ l enzyme in
50 μ l rxn.


8-16-96

* minis of suspected pEGFP-N2 ; minis on 1 ml culture ;
after add'n of 900 μ l 100% EtOH put in EtOH bath
for 15' and centrifuged for 10' w/ subsequent
aspiration ; resuspended in 30 μ l 1X TE

* susp. pEGFP-N2 digest w/ EcoRI

5 μ l DNA (susp. pEGFP-N2 ; pEGFP-N1 control ⁺ used 1 μ l)		
0.1 μ l RNase A	} 35x	3.5 μ l RNase A
1 μ l EcoRI		35 μ l EcoRI
1 μ l EcoRI Buffer		35 μ l EcoRI Buffer
2.9 μ l ddH ₂ O		101.5 μ l ddH ₂ O
10 μ l tot. vol.		175 μ l

* digest run on 1.5% gel :




EXPOSURE: 6/30 SEC. E-BOX 0 LACTE 255 GYFIVE 1.05
DATE: 08-16-1996 TIME: 18:44 ID# 569-01533 FILE: (UNTITLED)

* recultured #4, 6, 19, 22,
23 in

3ml Btamp

* Qiagen column on
#4, 6



#4
~180 ng
#6
~60 ng

SEC. E-BOX 0 LACTE 255 GYFIVE 1.05
TIME: 18:44 ID# 569-01533 FILE: (UNTITLED)

8-19-96

* pE6FP-N2 digested again w/ EcoRI & set up
dbl digestion w/ EcoRI & BspMI (will repprt.)

5 μ l DNA

1 μ l EcoRI (BspMI)

1 μ l EcoRI Buffer (buffer 2)

3 μ l ddH₂O

10 μ l tot. vol.

* digested pFC2 - 482 rglcZ w/ HindIII & BglII

30 μ l DNA

5 μ l Buffer B (Boehringer Buffer)

2.5 μ l HindIII

2.5 μ l BglII

10 μ l ddH₂O

50 μ l tot. vol.

all @ 37° 4 hr.

* second set dbl digestion on # 4, 6, 8, 12, 14, 15, 21-23

15 μ l DNA

1 μ l RNase A

2 μ l BspMI

10 μ l Buffer 2

72 μ l ddH₂O

100 μ l tot. vol.

incubate 37° overnight

10x

* digested pFC2 - 482 rglcZ run on 1% agarose
and purified ~730 bp fragment:

* resuspended in 50 μ l TE

8-20-96


* Digestion w/ Bsp MI & Eco RI to check for intron orientation:

* yesterday's Bsp MI rxn. phenol / chloroform extracted (1x)
ppt. w/ cold EtOH; resuspended in 10 μ l 1x TE;
digested w/ Eco RI as follows:

10 μ l DNA	
1 10 μ l Eco RI	} 15x
2 μ l buffer	
0.1 μ l RNase A	
6.9 μ l ddH ₂ O	
20 μ l tot. vol	
	15 μ l Eco RI
	20 μ l buffer
	1.5 μ l RNase A
	103.5 μ l ddH ₂ O

* digested 4 hr 37°

* run on 1.5% agarose gel; culture #4 & 23



* can see 905 bp band on #4 (tip) and can be discerned for #8, 22, 23; bad digestion w/ Bsp MI since can't see 250 bp band for negatives # 12, 14, 15, 21

EXPOSURE - 8.0 SEC. BLACK 0.75 WHITE 255.0. DATE: 08-20-1996 TIME: 12:54. FILE: (UNIT-23)

8-21-96

* digested pEGFP-N2 w/ Hind III & Bgl II so com
 ligate pFC2 -482 rglc Z :
 10 μ l DNA (~10 μ g)
 2.5 μ l Hind III
 2.5 μ l Bgl II
 5 μ l Buffer B
 30 μ l ddH₂O
 50 μ l tot. vol. @ 37° for 4 hr.



* run on 0.8% agarose gel ;
 Qiagen gel extract vector & resuspend
 in 50 μ l

* setup ligations w/ -482 rglc Z 73C
 fragments

cutEGFP	insert	ligase	lig. buffer	ddH ₂ O
1 μ l	2 μ l	1 μ l	1 μ l	5 μ l
↓	5 μ l	↓	1.5 μ l	6.5 μ l

* control = 1 μ l cutEGFP + 9 μ l ddH₂O

* ligations incubate @ 15° overnight

8-22-96

* plated ligation rxns onto amp + LB plates (100 μ l ea.)

8-23-96

* started 24 cultures in ~3 ml LB + amp

* picked from "2 μ l insert" plate

* did not observe in transformation efficiency

.10/1. mmis on cmv. 40 kip. rg/c. w GFP. V W RE (Janet)




Lane 7

FD 2315

Page 12

Signal: G:62 A:72 T:39 C:22

ACTTCCCTNAGAGCNCATACNCAATTTTNNFNDCTNTTCCCTT
340 360 380



Sample 08
Dye Terminator (AnyPrimer)
Lane 4

Signal: G:62 A:71 T:37 C:22

Points 60 6184 Base 1: 604
Instrument #907443
FD 2318

Wed, Aug 21, 1996 2:59 PM Page 12
X: 0 to 6947 Y: 0 to 1200
Spacing: 10.51

TTTCTCANNCCATTANGT 20 NCCGAGGCGACTCAANICTATCCTGCACTNTCCANG 80 ATNCCAGTTATCCANCJCCGATCCGGACATGTTNCCANAGGQAGGANCCT 120 CACCGGAGGTACC

CANGCANNNTGANGTGGACNCCNAGNTAAAGGGGACNAGTTCAJCGATGCCGGCAATGGCA 280 CCGATNCCATGCTNCCAN 220 CTGANCCTGAAGTTCATCJGCA CACCGGNCAN 280

TGCCGCTGCCCTGGGCAACGCTCGAGAC AACC TGA CCACTACAGTGGCANTGGGCTCTGCCGCTAGCCCGACCA 320 CCACTAACCAACACANAGANATGNNNATTCGGCCCATGTCGGA 380

ANGGTTNTTCCATTNA GCGCCAC 400 CTTCTTCATN GAGANNAA 420 GCGCCCGCTAGCAA CAAACGCGCCGCGCAAGCTGAANNTGGAAGGAGGAGAAACCCCT 480 TAAATCGGCTCTAA 540

ANGGTTNTTCCATTNA GCGCCAC 400 CTTCTTCATN GAGANNAA 420 GCGCCCGCTAGCAA CAAACGCGCCGCGCAAGCTGAANNTGGAAGGAGGAGAAACCCCT 480 TAAATCGGCTCTAA 540



Model 373A
Version 1.2.0


Sample 08.1
Dye Terminator(AnyPrimer)
Lane 8

Points 697, 6184 Base 1: 692
Instrument: 7907443
FD 2316

Wed, Aug 21, 1996 2:59 PM Page 1 of 2
X: 0 to 6561 Y: 0 to 1200
Spacing: 10.93

Signal: G:120 A:125 T:73 C:35

TCATGATGCGCAAGTCGAATTAACTCACTAAAGGAAACAAAGCGGAGCTGGCTTTAAAAAGCCCTCCACACCTCCCTGAACTGAAAGGTAAATGAATGCAATTGTTTAACTTGGTTA





Model 373A
Version 1.2.0

Sample 08.1
Dye Terminator(AnyPrimer)
Lane 8


Signal: G:120 A:125 T:73 C:35

Points 692 6184
Instrument: 907443
FD 2316

Base 1: 692

Wed, Aug 21, 1996 2:59 PM
X: 0 to 6561 Y: 0 to 1200
Spacing: 10.93


acgagatcctttaaattccccctnagggcncgcgggcccctttgttgcg



Signal: G:118 A:125 T:73 C:34

MC ~~2817~~ FD 2316

TCATGATGATGACGCGAAATTAACGCGCACTAAAGCGACAAACGCGGAGCGCGCGC | GAGC | AAAA
10 20 30 40 50 60 70 80 90 100 110 120 130



Model 373A

Version 1.2.0

Sample 19

Dye Terminator{AnyPrimer}

Lane 19 &

Signal: G:352 A:485 T:305 C:158

Points 68 6184 Base 1: 69T

instrument, #907443

~~Control PGEM~~ FD 2326


Wed, Aug 21, 1996 2:59 PM Page 112

X: 0 to 6503 Y: 0 to 1200

Spacing: 11.03

$$\frac{0}{\wedge}$$

PCORT-100 journal





Model 373A

Version 1.2.0

BT AACCCCAAGACATNAANGANCACAAITGATTCATTCT

Sample 19
Dye Terminator(AnyPrimer)
Lane 19


Signal: G:352 A:485 T:305 C:158


Points 69
Instrument: 907443
Control pgEM

Base 1: 691

Wed, Aug 21, 1996 2:59 PM
X: 0 to 6503 Y: 0 to 1200
Spacing: 11:03

Page 1/2







Model 373A


Version 1.2.0

480 586 510 528
TTCACACCGC TAC CCC ACC AGA TNAANCAAG ACAATTCTTGGCANTCT

Sample 18.1
Dye Terminator (Any Primer)
Lane 18
Signal: G:358 A:486 T:319 C:157

Points 60 6184 Base 1: 690
Instrument #907443
FD 2326

Wed, Aug 21, 1996 2:59 PM Page 2 of 2
X: 0 to 6631 Y: 0 to 1200
Spacing: 10.82



8-15-90

* Maize's frag. → used 50ng vector
~~used 1.0 insert~~ 3.0 insert

* 482 rgs promoter (rat glucagon)

* pEGFP-N2 = do double digest only w/ pos.

carrier DNA: digest ; when ppt. w/ EtOH add 1.0 μl
tRNA then continue

* plasmid prep → dbe CsCl (2x) of pEGFP-N2

* also reculture & do 20- tip of pEGFP-N2

* always leave 1 μl of DNA from prep.

* transform 1 μl DNA for a plate

~~for~~

* for BamHI / HindIII digest use Buffer B ;
can phenol / chloroform purify or gel purify

* intron - inc. exp. level & thus sens. of system


*


8-24-96

* minus pEGFP-N2 -482 rglc 2-element w/ EcoRI
overnight

pick minus for pEGFPN2-482+rglc 2 element

9-25-9





EXPOSURE: 12/30 SEI. BLACK: 0 WHITE: 255 GAMMA: 1.05
DATE: 08-25-1998 TIME: 09:56 ID: 569-01698 FILE: (UNTITLED)